architect

EN	CoCr CORONARY STENT SYSTEM. INSTRUCTIONS FOR USE	3
ES	SISTEMA DE STENT CORONARIO CoCr. INSTRUCCIONES DE USO	9
PT	SISTEMA DE STENT CORONÁRIO CoCr. INSTRUÇÕES DE USO	15
IT	SISTEMA DI STENT CORONARICO CoCr. ISTRUZIONI PER L'USO	21
FR	SYSTÈME D'ENDOPROTHÈSE CORONAIRE CoCr. INSTRUCTIONS D'UTILISATION	27
DE	KORONARSTENT SYSTEM CoCr. GEBRAUCHSHINWEISE	33
RU	СИСТЕМА КОРОНАРНОГО СТЕНТА ИЗ КОБАЛЬТ-ХРОМОВОГО СПЛАВА (CoCr). ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ	39
CS	SYSTÉM KORONÁRNÍHO STENTU CoCr. NÁVOD K POUŽITÍ	45
EL	ΣΥΣΤΗΜΑ ΣΤΕΝΤ ΣΤΕΦΑΝΙΑΙΑΣ CoCr. ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	51
PL	SYSTEM STENTU WIĘNCOWEGO CoCr. INSTRUKCJA OBSŁUG	57

INSTRUCTIONS FOR USE
CoCr CORONARY STENT SYSTEM

1. Description of the Product

The LVD Biotech SL coronary stent is a stent made of a Chromium cobalt alloy named L605. The stent is supplied pre-mounted on a delivery system which will expand it, by means of the distal balloon, in the coronary artery to be treated.

The stent is intended to fit with different artery diameters by adjusting its open-cell design with alternated connecting bridges. The stent is manufactured means of laser cuttings in metallic tubes. Afterwards, it is submitted to several treatments in order to give a smooth and sleek finishing to the surface.

Scheme of the delivery system with the stent pre-mounted onto the balloon. The product is supplied with different diameter and length configurations.

The stent delivery system consists on a rapid-exchange balloon catheter. The catheter comprises an inflatable balloon at the distal part. The balloon has different diameters and lengths to fit with the different stent configurations. When the balloon is inflated, it dilates the stent and deploys it against the artery. Afterwards, the balloon is deflated and removed, and the stent remains permanently implanted.

There are marks on the catheter shaft to help the user to calculate the position of the catheter as it is advanced through the guiding catheter without fluoroscopy, so that, when the last mark disappears, the catheter is already near to the tip of the guiding catheter and about to enter the artery. The marker located nearest to the catheter adapter is for femoral guiding catheters and the one furthest away is for brachial guiding catheters.

The distal part of the catheter is covered in a durable hydrophilic coating to lubricate the catheter so that it can navigate through arteries.

The maximum guide wire passage diameter must not exceed 0.36mm = 0.014".

Useful length of the catheter is 142cm, while the total length is 150cm.

The stent is supplied in the following stent lengths and diameters:

		STENT DIAMETERS (mm)								
		2	2.25	2.50	2.75	3.00	3.00	3.50	4.00	4.50
STENT LENGTH (mm)	9	X	X	X	X	X	X	X	X	
	14	X	X	X	X	X	X	X	X	X
	19	X	X	X	X	X	X	X	X	X
	24	X	X	X	X	X	X	X	X	X
	29	X	X	X	X	X	X	X	X	X
	34	X	X	X	X	X	X	X	X	X
	39	X	X	X	X	X	X	X	X	X

To expand the stent, the proximal luer-lock hub is needed to be connected to an inflation device with a manometer. The balloon is inflated with the pressure shown by the manometer to expand the stent to predictable diameters. Two radiopaque markers are located at each end of the balloon in order to mark its length and help the user to know its position while it is inside the patient.

Contents

- One coronary stent system comprising the stent pre-mounted on its delivery system (balloon catheter). The stent is protected by a sheath and a stylet inserted at the guide insertion lumen, and it is inserted into a dispenser and in a sterile bag.
- One card with the compliance curve showing the working pressure range

INSTRUCTIONS FOR USE CoCr CORONARY STENT SYSTEM

- One implant card
- One leaflet with instructions for use

2. Indications

The device is indicated for increasing the internal diameter of an artery with the aim of improving blood flow in the following cases:

- Patients with symptomatic ischemic heart disease due to "de novo" stenotic and re-stenotic lesions located in arteries with diameters from 2 mm to 4.5 mm.
- Patients with ischemic disease at coronary bypass grafts, including saphene vein grafts.
- Patients with occlusive disease due to acute myocardial infarction located in arteries with diameters from 2 mm to 4.5 mm.

3. Contraindications

Antiplatelet and/or anticoagulation therapy is contraindicated.

Patients with clear diagnosed heavy metal allergy.

Patients with lesions that prevent complete angioplasty balloon inflation or proper stenting.

4. Warnings

- Stent implantation must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or, in its absence, in hospitals with easy access to other hospitals where this type of operation can be performed.
- The product should only be used by experienced physicians in hemodynamic and familiar with the percutaneous coronary intervention and stent implantation.
- The device is designed for single use in a single patient. DO NOT RE-STERILIZE OR RE-USE IT. REPROCESSING IS FORBIDDEN. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check its expiry date and do not use products which are past this date.
- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., according to the protocol of coronary intravascular catheter insertion.
- Use aseptic techniques when the product is taken out of its primary packaging.
- Do not dry with gauze.
- Do not expose the delivery device to organic solvents.
- Do not use oily or high viscous contrast media unsuitable for intravascular use.
- Do not use air or gaseous media to inflate the balloon, it must be inflated with a mixture of saline solution and contrast liquid (preferably 1:1).
- Select the suitable product size, in terms of diameter and length, according to the size of the lesion fluoroscopically observed.
- Advance the product on the guide wire by means of fluoroscopy imaging. Do not allow product advance without the guide wire inside it.
- Do not overlap stents with different compositions.
- Do not expose the product to hospital waste material.

5. Precautions

5.1. General precautions

- It is recommended to be stored in a cool, dry place (0°-40°C), away from direct sunlight.
- Inspect the packaging. If it is damaged, do not use the product.
- After its use, this product may be a biological hazard. Handle and dispose it of according to the accepted medical practices, and the pertinent local, state or federal laws and standards.

5.2. Precautions during product preparation

- Inspect the product before using it to ensure it is in perfect conditions, and discard any device not showing proper confidence.

INSTRUCTIONS FOR USE
CoCr CORONARY STENT SYSTEM

- Do not prepare or pre-inflate the delivery system before starting the procedure or deploying the stent at the lesion.
- Do not handle the stent with the fingers, due to this action may dislodge the stent from the balloon.
- Do not attempt to separate the stent from its delivery system; if any stent movement or unsuitable stent place or crimping is noticed, do not use it.
- Do not attempt to straighten a kinked hypotube (proximal part of the delivery system), it may get broken.

5.3. Precautions during product advance-positioning-stenting-withdrawal

- In case of multiple lesions treatment, begin with the most distal lesions, followed by the proximal ones.
- Always advance the device with fluoroscopic imaging. If any resistance is noted, do not continue until determining its cause. If the stent is not able to cross or reach the lesion, remove the whole system (including the guiding catheter and the guide wire) as a single unit.
- If any resistance is noticed during advance, stop the advance and determine its cause before continuing.
- Once stent expansion is started, do not intent to remove it or change its position.
- Do not manipulate, advance or retract, either the catheter or the guide wire, while the balloon is inflated.
- Do not exceed the maximum recommended pressure indicated on the label and the attached compliance curve.
- If any resistance is noted during lesion access, delivery system withdrawal or withdrawal of a stent which could not be implanted, the whole system should be removed as a single unit:
 - Do not withdraw the delivery system into the guiding catheter
 - Place the proximal balloon marker distally to the tip of the guiding catheter
 - Advance the guide wire towards into the coronary anatomy as far as safely possible
 - Tighten the haemostatic valve so it properly grabs the guiding catheter and the balloon
 - Remove the guiding catheter and the balloon catheter together, as a single unit
 - Finally, remove the guide wire, or in case an implantation is wished to begin again, maintain the guide wire and placed the rest of the devices.
- Great care must be exercised immediately after implantation and if it was needed to re-cross it, not to push the stent with guide wires or balloons not to destructure it.

5.4. Precautions: MRI¹ Safety:

Non-clinical trials have demonstrated safety on obtaining L605 type chromium cobalt stent images by magnetic resonance (MR) under certain conditions (MR conditional). MRI test conditions used to evaluate these devices have been: for magnetic field interactions, a static magnetic field strength of 3.0 Tesla, with a maximum spatial gradient magnetic field of 720 gauss/cm; for MRI-related heating, a maximum whole body averaged specific absorption rate (SAR) of 1.54W/kg for durations of 15 minutes within a 1.5 Tesla field and an average specific absorption rate (SAR) of 2.8W/kg for a duration of 15 minutes within a 3.0-Tesla field. Temperature increase for a single stent was no greater than 0.8°C, while for an overlapped pair it was no greater than 1.0°C.

No clinical or non-clinical trials have been conducted to rule out the possibility of stent migration at an over 3.0 Tesla strength field. It is not recommended to perform MR tests with an over 3.0 Tesla strength until 7 days after implantation.

MR imaging quality may be affected if the area of interest is adjacent or very close to the stent place.

6. Possible Adverse Effects / Complications

Possible associated adverse effects and/or complications that may arise before, during or after the procedure include the following:

- Death
- Acute myocardial infarction
- Unstable angina

¹ These data were obtained bibliographically, by searching for marketed stents with the same composition.

INSTRUCTIONS FOR USE
CoCr CORONARY STENT SYSTEM

- Ictus / embolism / thrombosis
- Emergency bypass surgery
- Endocarditis
- Allergic reaction to the metals or contrast media
- Coronary artery spasm
- Hemodynamic deterioration
- Arrhythmia
- Aneurysm or pseudo-aneurysm
- Infections
- Total occlusion of the artery / thrombosis / stent embolization
- Reocclusion of the treated area: re-stenosis
- Perforation or dissection of the treated area
- Arteriovenous fistula
- Hypotension/hypertension
- Local haemorrhaging with haematoma at the insertion site

7. How to Use

7.1. Required equipment

- Heparinised normal saline solution
- Contrast medium (use contrast media indicated for intravascular use)
- Proper medical prescription to the procedure
- 0.014" guide wire: do not use any other size of guide wire
- Introducer with haemostatic valve whose size is equal or larger than 5F (except for 4 and 4.5 mm stent diameters, where 6F is recommended). Do not use smaller sizes as these might damage or leave useless the catheter
- Guiding catheter
- Haemostatic valve
- Three-way stopcock
- Extension lines
- Inflation device (manual pump with built-in manometer)
- Several standard 10-20cc syringes with saline solution for washing the system.

7.2. System Preparation

- The stent diameter must coincide with the reference vessel diameter, as its length with lesion length. It is preferably to be a slightly longer to ensure full lesion coating.
- Remove the catheter from the protective dispenser. Check that it is the right size. Remove the protective pod and stylet from the balloon. Wash the guide lumen with sterile saline solution by means of injection with a syringe from the tip until it outcomes from the guide wire insertion port. Do not dry with gauze.
- Inspect the stent to ensure it is located between the radiopaque marks; if it is noted that it has moved, do not use it.
- Submerge the catheter in sterile saline solution to activate the hydrophilic coating.
- Attach a three-way stopcock to an extension line and the extension line to the luer hub of the delivery system. Optionally, the inflation pump can be connected to the three-way stopcock.
- Close the air passage through the balloon.
- Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
- Open the three-way stopcock connection between the syringe and the catheter.
- Withdraw the syringe plunger upwards allowing the air bubbles exiting to the liquid (it is recommended that the syringe be in a vertical position).
- When bubbles stop coming in, leave a meniscus of solution in the key and close the three-way stopcock at the catheter part. Remove the syringe.
- Connect the inflation device with 1/3 of 1:1 contrast medium and saline solution mixture (duly

INSTRUCTIONS FOR USE CoCr CORONARY STENT SYSTEM

purged), to the three-way stopcock and the extension line. Do not apply pressure as this might dislodge the stent.

DANGER: If bubbles do not stop coming from the catheter into the syringe after 1 minute of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly sealed. If bubbles continue to be seen after checking the connections, do not use the device. Return the device to the manufacturer or distributor for its inspection.

7.3. Implantation Technique

Preliminary steps

- Insert the introducer according to the manufacturer's instructions.
- Connect the haemostatic valve to the guiding catheter hub, maintain the valve closed. Insert the guiding catheter through the introducer according to the manufacturer's instructions and make sure it is located at the coronary artery entrance. Inject contrast liquid through the guiding catheter to ensure it is correctly positioned giving access to the coronary artery.
- Open the haemostatic valve for 0.014" guide wire passage through. Once the guide crosses the key valve, seal it to avoid blood loss. Place the 0.014" guide wire through the lesion according to percutaneous intervention techniques, by means of fluoroscopy techniques to determine its position all times. Use a contrast liquid infusion to determine that the guide wire is properly positioned.

Implantation technique

When the lesion diameter is so small, calcified or as tortuous that may compromise stent access, its pre-dilatation with a coronary dilatation balloon catheter should be considered. Follow the instructions of the balloon catheter manufacturer for lesion dilatation. Then, proceed to stent implantation.

- Insert the proximal end of the guide wire into the distal tip of the delivery system.
- Open the haemostatic valve again to allow stent passage through the valve. Once the stent has been run through it, close the valve. If any resistance is noted, do not advance the delivery system through the haemostatic key. Avoid the haemostatic key valve damaging the catheter shaft; which may subsequently affect the balloon in-deflation.
- Carefully advance the stent system on the guide wire through the guiding catheter until it reaches the lesion to treat, by means of fluoroscopic imaging. Ensure that the guide wire comes out of the balloon catheter guide insertion port (approximately 25cm from the tip of the balloon catheter).
- Place the stent across the lesion situating the radiopaque markers on both sides of the lesion. Use high-resolution fluoroscopy to verify that the stent has not been damaged or shifted during its advance before its expansion.
- Apply pressure for 15-30 seconds depending on the lesion, to completely expand the stent. Do not exceed the rated burst pressure (RBP) indicated at the label and the compliance curve.
- For an optimal expansion, the stent must be in contact with the artery, otherwise the stent could move when the balloon is removed. The balloon diameter must be slightly higher than the vessel diameter to allow stent elastic recoil once the balloon is deflated. Consult the compliance table for proper expansion, due to too high size may also induce artery dissection.
- Withdraw the plunger of the inflation device to deflate the balloon. Keep up the negative pressure for between 15 and 30 seconds depending on the balloon size. Ensure that the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire still placed, withdraw the catheter from the lesion. Keep the guide across the dilated stenosis.
- Carry out an angiography through the guiding catheter to confirm suitable stenting.
- Keep the guide wire through the dilated stenosis for 30 minutes after angioplasty. In case of any thrombus observed around the stent, administer thrombolytic agents. When the angiography confirms dilatation, carefully withdraw the guide wire and the balloon catheter through the connector.
- Withdraw the guiding catheter through the introducer

INSTRUCTIONS FOR USE
CoCr CORONARY STENT SYSTEM

- Leave the introducer placed until the haemodynamic profile returns to normality. Perform suture as usually.
- Further dilatation of stented segments

If the stent size results inadequate with respect to vessel diameter, stent post-dilatation technique may be applied with a larger balloon. To do this, a balloon catheter is inserted into the stented area (it is not needed to be compliant) and dilatation is performed according to the manufacturer's recommendations.

As a further recommendation, a "small stent" must not be post-dilated over 3.30mm, a "medium stent" must not be post-dilated over 4.25mm and a "large stent" must not be post-dilated over 5.25mm.

Withdrawing an unexpanded stent

If the stent does not cross the lesion, it may need to be removed without expanding. If any resistance is not detected, check that the guiding catheter is coaxially situated with respect to the stent and carefully withdraw the stent towards the guiding catheter. If any unusual resistance is noticed during its withdrawal, extract the whole system as a single unit by means of fluoroscopic imaging according to the indicated instructions in part 5.3 of this leaflet.

8. Warranty

The product and all its components are designed, manufactured, tested and packed taking all reasonable precautions. LVD Biotech guarantees the product up to its expiry date providing the packaging is not broken, tampered with or damaged.

STANDARDISED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES
(EN 980-2008 STANDARD)

			
Catalogue reference	Batch no.	Expiry date	Keep away from sunlight
			
Manufacturer	Sterilized with ethylene oxide	Limitation in storage temperature	Keep dry
			
Consult instructions for use	Do not re-use	Do not re-sterilize	Do not use if the pack is damaged

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© iVascular - cod. MP45005 rev2 Jun 2013